

Appendix 5. Annotated list of a new radiolaria- and diatom-slide preparation

Twenty-seven sub-samples from H.M.S. Challenger Collection deposited in the Natural History Museum in London were provided with the project. For the future studies on microfossils in the *Challenger Expedition* we made slides from the samples with the following procedure, and these should be deposited in the five institutions: NHM in London, MN in Berlin, Tohoku Univ., Utsunomiya Univ. and NMNS in Tokyo.

Additional Nineteen sub-samples from H.M.S. Challenger Collection housed in the NHM, London were provided to the project. Regarding these samples, only radiolarian slides were made and a 19 slide set was distributed and deposited in the same five institutions.

Methods for Preparation of Radiolarian Slides

- Disaggregation** The sample was placed in a 300ml beaker with the boiled solution of hydrogen peroxide (H₂O₂, 10%, 150ml) and hydrochloric acid (HCl, 5%, 30ml).
- Washing** After boiling for about 40-60 minutes, a suspension was added by water and poured into a 45µm-sieve and washed.
- Cleaning** Washed residue was placed in a 300ml beaker with the warmed (60-80°C) solution of sodium hexametaphosphate (1-2%, 200ml) to disperse clay grains within a shell for 40 to 60 minutes.
- Drying** The cleaned residue was washed upon a 45µm-sieve and dried.

Preparation of Strewn Slides

- Dispersion** The dried residue was placed on a glass slide to which was already applied a thin glue of a gum tragacanth. Then, the residue was moisted by a wet air and then fixed on the slide glass after it dried. Turning over the slide glass and an excess residue was freed.
- Mounting** A single drop of xylen was poured on the center of the residue and then was mounted on cover glass with a Canada Balsam.

Methods for Preparation of Diatom Slides

Chemical & Physical Cleaning

- Drying** sediment samples dried at 60°C for 24 hours
- Weighing** about 0.5 - 1g were weighed at an accuracy of 0.001g
- Disaggregation** the sample was placed in a 200ml beaker with the boiled solution of hydrogen peroxide (H₂O₂, 15%, 20ml) and hydrochloric acid (HCl, 5%, 1-5ml)
- Centrifugation** after boiling for about 20 minutes, a suspension was poured into a centrifuge tube and filled with water
the suspension was centrifuged for 2 minutes at a speed of 1200rpm
fine material in suspension was carefully removed by decantation, and the tube was filled again with water
this procedure was repeated 5 times
- Suspension** the residue was diluted to 10ml with water and stored in a glass vial

Preparation of Slides

- Dilution** the suspension was placed in a short test tube by the aid of an automatic micropipette, and diluted by water to obtain a suspension of proper density of diatom valves
- Drying** the diluted suspension was placed on a square cover glass, 18 x 18mm, and dried at 45°C
- Mounting** the sample was mounted on a slide glass with StyraX